


# Dissolved Oxygen (DO)

- ❑ Dissolved oxygen (DO) is the *amount of oxygen dissolved* in a given quantity of water at a particular temperature and atmospheric pressure
- Dissolved oxygen gets into the water by diffusion from the atmosphere. Aeration of the water as  $O_2$  is waste product of photosynthesis
- It is an important parameter in assessing water quality because of its influence on the organisms living within a body of water
- Amount of oxygen dissolved is expressed as mg/L or ppm




- **Biochemical Oxygen Demand (BOD):** Amount of dissolved oxygen required by bacteria/other microorganisms to break down organic material present in a given water sample at certain temperature over a specific time period.
  - To determine BOD, a measured volume of waste water is taken and its DO is measured. It is then incubated for 5 days at 20 °C and its DO is determined again.
- **Chemical Oxygen Demand (COD):** An indicative measure of the amount of oxygen that can be consumed by reactions in a measured solution.
  - It is commonly expressed in mass of oxygen consumed over volume of solution which is milligrams per litre (mg/L).
  - A COD test can be used to easily quantify the amount of organics in water.

## Causes of Low Dissolved Oxygen

- Low dissolved oxygen (DO) primarily results from **excessive algae growth caused by Nitrogen/phosphorus** which generally originate from discharges from municipal and private wastewater treatment, cropland and urban storm water runoff, and natural decay of vegetation
- After algae die, decompose process consumes dissolved oxygen. This can result in insufficient amounts of dissolved oxygen available for aquatic species.
- Die-off and decomposition of submerged plants also contributes to low dissolved oxygen.






- Aeration
- Water temperature
- Organic wastes
- Aquatic plant populations
- Photosynthetic activity of the water
- Respiration of animals and plants
- Speed of water flow
- Altitude (atmospheric pressure)
- Human activities
- Salt concentration: The solubility of oxygen decreases with increase in concentration of the salt.
- Other dissolved gases concentrations



## # Estimation of DO by Winkler's method

- DO is important with respect to the species of aquatic life, It is also measure of its ability to oxidize organic impurities in water
- DO test is used to control the amount of oxygen in boiler feed water
- DO test is helps to assess the raw water quality and to keep a check on stream pollution
- It is one of the most useful titrations involving iodine
- In order to avoid loss of oxygen from the water sample, it is fixed by  $\text{MnSO}_4$

- △ The principle involved in this methods of determination of DO is to bring about the **oxidation of potassium iodide (KI) to iodine with dissolved oxygen** present in water sample after adding  $\text{MnSO}_4$ , KI and KOH
- △ The reaction with **manganese(II) hydroxide which is converted rapidly and quantitatively to manganese(III) hydroxide**. Here  $\text{MnSO}_4$  acts oxygen carrier to enable the dissolved oxygen in molecular form to take part into the reaction
- △ On **acidification**, the **manganese reverts back to the divalent state** and an **equivalent amount of iodine is liberated form the KI** present
- △ The liberated  $\text{I}_2$  is **titrated against standard sodium thiosulfate (hypo) solution**, using starch as indicator
- △ This means that **4 mol thiosulphate** is equivalent to **1 mol dissolved oxygen**



---

## □ Procedure:

- **Standardization of  $\text{Na}_2\text{S}_2\text{O}_3$**
- **Estimation of DO using standardized  $\text{Na}_2\text{S}_2\text{O}_3$**



# DO estimation by Winkler's method

## □ Procedure:

### ➤ Standardization of $\text{Na}_2\text{S}_2\text{O}_3$

- The secondary standard solution of sodium thiosulphate is standardized by titrating with a primary standard potassium dichromate using **starch** as **indicator**
- Color change occurs from **straw yellow** to **blue** to **colorless**

