Calculate Total hardness of the given hard water sample with following data:
 [i] 50 mL of standard hard water sample consume 42 mL of EDTA solution.
 [ii] 50 mL of hard water sample consume 20 mL of EDTA solution before boiling.
 Standard hard water is prepared by dissolving 1.0 gram of CaCO₃ per Litre.

1000 mL of Standard hard water contains 1 g of CaCO₃

⇒ 50 ml of standard hard water contains = 50/1000 g of CaCO₃

Given that, 42 mL of EDTA is consumed by 50 mL of standard hard water

42 mL of EDTA = 50/1000 g of CaCO₃

∴1 mL of EDTA = [50/(42*1000)] g of CaCO₃

⇒ EDTA consumed by sample hard water = 20 mL
 So, 20 mL of EDTA = [(50*20)/(42*1000)] g of CaCO₃
 Hence, 50 mL of sample hard water contains [(50*20)/(42*1000)] g of CaCO₃
 ∴ 1000 mL of sample hard water = [(50*20*1000)/(42*1000*50)] g of CaCO₃ = 0.476 g of CaCO₃
 i.e. Total hardness of sample hard water = 0.476 g = 476 mg of CaCO₃ /liter = 476 ppm

Problem 2 for DA



- Calculate Total, Temporary and Permanent hardness of the given hard water sample with following data:
 - [i] 50 mL of standard hard water sample consume 42 mL of EDTA solution.
 - [ii] 50 mL of hard water sample consume 19.9 mL of EDTA solution before boiling.
 - [iii] 50 mL of hard water sample consume 8.9 mL of EDTA solution after boiling.
 - Standard hard water is prepared by dissolving 1.5 gram of CaCO₃ per Litre.